



POLITÉCNICA

“Ingeniamos el futuro”

INSE

Práctica: Smart Health

Reconocimiento de actividades físicas basado en acelerómetro 3D

Cristina Luna Jiménez,

Fernando Fernández-Martínez,

Jose Manuel Pardo Muñoz,

fernando.fernandezm@upm.es

Contents

Introduction


HW: Smartphone

SW architecture

Assignment

Weka Experimenter





Embedded Intelligence for Health Care and Wellbeing



a



b



c



d



e



f

Sensors embedded
in everyday
surroundings

a) armband

b) ECG monitoring on a
car seat

c) neonatal monitoring

d) ear worn activity
recognition

e) Nike+iPod for
monitoring running

f) wheelchair motion
and paralympic
athlete performance

Medical Internet of Things (MIoT)

- Integration of medical devices in a network connection
- Network can be managed from the web
- Provide information in real time
- Communication: person to person (P2P) & machine to machine (M2M)
- Allow interaction between health professionals & patients
- MIoT can be seen from 3 paradigms:
 - Internet-oriented middleware
 - Things sensors oriented
 - Knowledge-oriented semantics

Embedded Intelligence in MIoTT

- Growth at a high rate exceeding 7%
- Estimated Revenues by 2020 \$2.2 trillion
- Healthcare is one of the leading industries
 - Source: 3rd WHO Global Forum on Medical Devices ([link](#))


What is it?

- **Embedded intelligence** is the ability of a product, process or service to monitor its
 - Operational performance,
 - usage load,
 - environment
- Goals:
 - enhance performance and lifetime,
 - increase quality and
 - improve customer satisfaction

How does it help?


- Monitor health and usage of products to extend their performance and lifetime
- Improve market appeal and acceptance of products
- The ability for a service, system or product to be used by ageing and people with special needs
- Address skills shortages in limited resources
- Enabling new revenue opportunities

Example: Architecture of System



Health monitoring: Physical activity recognition

- System architecture:



Health monitoring: Demo



Contents

Introduction

HW: Smartphone


SW architecture

Assignment

Weka Experimenter




Sense HAT




- Add-on board for Raspberry Pi
 - Made especially for the Astro Pi mission
- 8×8 RGB LED matrix
- Five-button joystick
- Sensors:
 - Gyroscope
 - **Accelerometer**
 - Magnetometer
 - Temperature
 - Barometric pressure
 - Humidity
- Python library

Smartphone



- Sensors:
 - Gyroscope
 - **Accelerometer**
 - Magnetometer
 - Camera
 - ...



Contents

Introduction

HW: Smartphone

SW architecture

Assignment


Weka Experimenter



Software tools

Colab allows you to write and execute Python in your browser with:

- No configuration
- Free access to GPUs
- Easy sharing



Drive: file storage and synchronization service developed by Google



Python: is an interpreted, high-level and general-purpose programming language

- 1) Download [this](#)
- 2) Upload it to your personal Google Drive account
- 3) Double-click on it to launch Colab

Introduction to Colab

Setting GPU environment

<https://colab.research.google.com/>



Introduction to Colab

The screenshot shows the Google Colab interface. On the left, the 'Files' pane displays a directory structure including 'bin', 'boot', 'content', 'logs', 'models', 'sample_data', 'samples_app', 'adc.json', 'df_test.csv', 'df_test_nominal.csv', 'df_train.csv', 'df_train_nominal.csv', 'octave.zip', 'datalab', 'dev', 'etc', 'home', and 'lib'. A green box highlights the 'Files' icon in the top left of this pane, with an arrow pointing to a text box labeled 'List of folders/files'. The main area contains a code editor with two cells. The first cell contains a list of package installation commands. The second cell contains Python code for installing Octave and using PyDrive. A green box highlights the 'Run' button (a play icon) in the first cell, with an arrow pointing to a text box labeled 'Run a code cell'. Below the code editor, a green box highlights the '+ Add new cell' button, with an arrow pointing to a text box labeled 'Create new cell (Text or Code Cell)'. A tooltip for this button shows 'Añadir celda de código Ctrl+M B'.

List of folders/files

Run a code cell

Create new cell
(Text or Code Cell)

Software architecture

- Main colab script:
 - Practice SmartHealth
- Project scenarios:
 - Raw data + Decision Tree
 - Feature Extraction + Decision Tree
 - Raw data + Temporal Model

Índice



Practice: Smart Health

Introduction



Environment Set-up



Functions

Weka Training-Validation-Testing

Feature Extraction

LSTM Functions

Common

Baseline

Download datasets

Raw data + Decision Tree Scenario

MFCC Feature Extraction + Decision Tree Scenario

Raw data + Temporal Model (LSTM)

Systems Comparison

Assignment

Appendix

Upload files (or datasets)

Download files (datasets, models...)

Save files in your Colab session

Save and Load networks

Dataset as ARFF

Physics Toolbox Suite

Gesture recognition

- Basic (table 😊) tennis techniques or strokes:



Drive (1)





Backhand (2)




Serve (3)


Lob (4)



Gesture recognition



- Based in a 3D accelerometer
 - **Windowed analysis** of acceleration values
 - Different window lengths
 - Overlapped windows
 - Low sampling frequency: 50 Hz
- Simple front-end:
 - Only 60 mfcc-plp features
 - 20 features per channel
 - 6 functionals (mean & stdev)
 - Classifiers: Random Forest, LSTM




Feature Extraction




- Scientific Programming Language
- Powerful mathematics-oriented syntax with built-in plotting and visualization tools
- Free software, runs on GNU/Linux, macOS, BSD, and Windows
- Drop-in compatible with many Matlab script

```
# Create an evenly-spaced vector from -10..10  
x = -10:0.1:10;  
y = sin(x); # y is also a vector  
plot(x, y);  
title("Simple 2-D Plot");  
xlabel("x");  
ylabel("sin (x)");
```



Feature Extraction (2)




Feature Extraction (3)



```
FeatureExtraction (AccX_array, AccY_array, AccZ_array (  
    200, # WINDOW_SIZE: K (4 seconds)  
    150, # OVERLAP: K - Q (3 seconds)  
    False)
```

Windowing is a means to **stationarize** signals.
Tools like **Fourier-based techniques** are well-suited to stationary signals.



Feature Extraction (4)

```
from oct2py import octave # Oct2Py: Python to GNU Octave Bridge
```

```
def FeatureExtraction (AccX_array, AccY_array, AccZ_array,
window_size, overlap, debug):
```

```
    aux_AccX_array = AccX_array[np.newaxis, :].T
```

```
    aux_AccY_array = AccY_array[np.newaxis, :].T
```

```
    aux_AccZ_array = AccZ_array[np.newaxis, :].T
```

```
    ...
```

```
    mfcc_plp = octave.calcula_features_mfcc_plp_online_new(aux_AccX
aux_AccY_array,          aux_AccZ_array, window_size, overlap, False)
```

```
%30plp+30mfcc
```

```
    functionals = [np.mean(aux_AccX_array), ..., np.std(aux_AccX_array), ...]
```

```
%30plp+30mfcc+6functionals=66 features)
```

```
/usr/local/lib/octave/calcula_features_mfcc_plp_online_new.m
res = np.append((np.asarray(mfcc_plp, ...)), (np.asarray(functionals, ...)))
```

```
function OUTPUT =
```

```
calcula_features_mfcc_plp_online_new (AccXorig, AccYorig, AccZorig,
win_size, overlap)
```

```
    ...
```

```
    OUTPUT = frastaplp3D (AccX', AccY', AccZ', win_size, step); %30 features
(10+10+10)
```

```
    OUTPUT = [OUTPUT fmfcc3D (AccX', AccY', AccZ', win_size, step)]; %30 features
(idem)
```

Oct2py
runs m-files
from python



Generated dataset



dataset_test_mfcc

...	X_PLP_param1	X_PLP_param2	X_PLP_param3	X_PLP_param4	X_PLP_param5	Z_MFCC_param60	X_mean_param	Y_mean_param	Z_mean_param	X_std_param	Y_std_param	Z_std_param	user	class
0	27.958436	1.792896	0.926759	0.278980	0.390426	0.974356	6.562834	6.638258	0.282533	2.511517	3.209461	3.624661	15	DRIVE
1	28.394158	1.744596	0.834124	0.150108	0.108082	-0.998596	5.054386	6.178057	-1.635454	5.180306	6.246923	5.153387	15	DRIVE
2	28.271179	1.223328	1.108851	0.755840	0.342047	-0.809730	4.535134	6.343756	-2.171127	5.488920	7.083468	5.056593	15	DRIVE
3	29.710130	1.286819	-0.219463	0.771217	-0.053580	-2.450957	4.061923	6.211137	-0.709630	5.387086	7.425837	6.012131	15	DRIVE
4	29.605985	1.300449	-0.039957	0.903575	0.167716	-1.657046	3.639123	6.371750	1.139097	5.470489	8.576052	7.393074	15	DRIVE
...
3835	28.423973	1.272255	0.619025	0.892203	-0.743771	-1.966657	-2.958946	4.502018	5.495929	4.909814	3.949686	7.236630	20	LOB
3836	28.460169	1.431241	0.651814	0.437103	-0.689868	-2.727839	-2.926883	4.677041	5.925444	4.919016	4.100518	6.692011	20	LOB
3837	28.448223	1.415562	1.002449	0.808588	-0.635030	-0.990033	-4.033475	5.615216	6.889087	4.730290	4.327907	6.212749	20	LOB
3838	29.302487	1.041828	0.690837	0.819883	-0.069165	0.479920	-4.806968	6.434846	7.822439	4.194530	4.210162	5.109674	20	LOB
3839	29.445724	0.845920	0.617501	0.816422	-0.001756	0.715593	-5.614426	7.062794	8.811072	3.497779	3.999955	3.995048	20	LOB





Create Model Weka

/Functions/Weka Training-Validation-Testing


```
from weka.classifiers import Classifier
def CreateClassifier(classname,options):
    model = Classifier(classname=classname, options=options)
    return model
```

/Feature Extraction + Decision Tree Scenario

options_RandomForest

```
=["-P","100", "-I", "100", "-num-slots","1", "-K", "0","-M","1.0", "-V", "0.001", "-S", "1"]
```

```
model_architecture_mfcc = CreateClassifier(
    classname="weka.classifiers.trees.RandomForest" ,
    options=options_RandomForest )
```



Create Model TF

/Functions/LSTM Functions

```
from weka.classifiers import Classifier

def CreateLSTM(neurons_first_layer, n_steps, dimension, loss, optimizer, metrics):
    model_lstm = Sequential()
    model_lstm.add(LSTM(neurons_first_layer, input_shape=(n_steps, dimension),
                        return_sequences=False))
    model_lstm.add(Dense(4, activation='softmax'))
    model_lstm.compile(loss=loss, optimizer=optimizer, metrics=metrics)
    return model_lstm
```

/Raw data + Temporal Model

```
model_lstm_ES = CreateLSTM(neurons_first_layer = 200, n_steps = 200, dimension =
3, loss='categorical_crossentropy', optimizer='rmsprop', metrics=['accuracy'])
```


Train Model Weka



/Functions/Weka Training-Validation-Testing

```
def TrainModel(train_dataset, classifier, output_path_model):  
    #Indicate where is class label  
    train_dataset.class_is_last()  
    # train classifier  
    classifier.build_classifier(train_dataset)  
    # save classifier object  
    print("\n[TrainedModel]\n\n")  
    serialization.write(output_path_model, classifier)  
    return classifier
```




```
#Train and save model  
out_path_model_mfcc = "models/trained_TreeRandomForest_mfcc.model"  
trained_model = TrainModel(dataset_train_mfcc, model_architecture_mfcc, out_path  
_model_mfcc)
```



Train Model TF

```
#Train model
```

```
model_lstm_ES.fit(X_train,y_train,batch_size=16,epochs=10,  
                  validation_data=(X_validation, y_validation),...)
```





Recognition Weka

/Functions/Weka Training-Validation-Testing

```
def Recognize(trained_classifier, test_dataset):  
    print("# - actual - predicted - error - class distribution" )  
    for index, inst in enumerate(test_dataset):  
        pred = trained_classifier.classify_instance(inst)  
        dist = trained_classifier.distribution_for_instance(inst)  
        print("%d - %s - %s - %s - %s" %  
              (index+1,  
               inst.get_string_value(inst.class_index),  
               inst.class_attribute.value(int(pred)),  
               "yes" if pred != inst.get_value(inst.class_index) else "no",  
               str(dist.tolist())))
```

#SEE PREDICTIONS:

```
Recognize(trained_model, dataset_test_weka)
```

Recognition TF

/Functions/Weka Training-Validation-Testing

```
def Recognize_LSTM(trained_LSTM, test_parameters, test_labels):
    y_pred = trained_LSTM.predict(test_parameters)
    print("# - actual - predicted - error - class distribution" )
    for index in range(0, len(test_labels)):
        pred = Convert2nominal_pingPong (np.argmax(y_pred[index], axis= -1)+1)
        true_val = Convert2nominal_pingPong (np.argmax(test_labels[index], axis= -1)+1)

        print("%d - %s - %s - %s - %s" %
              (index+1,
               true_val,
               pred,
               "yes" if pred != true_val else "no",
               str(y_pred[index])))
```

#SEE PREDICTIONS:

```
Recognize_LSTM(model_lstm_ES, X_test, y_test)
```

Evaluation Weka



```
/Functions/Weka Training-Validation-Testing  
from weka.classifiers import Evaluation  
def EvalClassifier (trained_classifier, train_dataset, test_dataset):  
    evl_train = Evaluation(train_dataset)  
    evl_test = Evaluation(train_dataset)  
    print("-----RESULTS IN TRAINING-----" )  
    evl_train.test_model(trained_classifier, train_dataset)  
    print(evl_train.summary())  
    print("-----RESULTS IN TEST-----" )  
    evl_test.test_model(trained_classifier, test_dataset)  
    print(evl_test.summary())
```

```
#EVAL TRAINED MODEL WITH MFCC FEATURES  
EvalClassifier(trained_model, dataset_train_mfcc, dataset_test_mfcc)
```




Evaluation TF

```
/Functions/Weka Training-Validation-Testing  
from weka.classifiers import Evaluation  
def EvalLSTM(y_labels, y_pred):  
    y_test_df = pd.DataFrame({ 0:y_labels[:, 0], 1:y_labels[:, 1], 2:y_labels[:, 2], 3:  
y_labels[:, 3]})  
    y_test_df_idx=y_test_df.idxmax(axis= 1)  
    print(confusion_matrix(y_test_df_idx, y_pred))  
    print(classification_report(y_test_df_idx, y_pred))  
    print('Accuracy %0.4f' % accuracy_score(y_test_df_idx, y_pred))
```

```
#Eval model on training/test sets  
print("---- TRAINING RESULTS: ----" )  
EvalLSTM(y_train, y_pred_train_ES)  
print("---- TEST RESULTS: ----" )  
EvalLSTM(y_test, y_pred_test_ES)
```


Systems Comparison



Contents

Introduction

HW: Smartphone

SW architecture


Assignment

Weka Experimenter




Assignment

- **STANDARD: Re-run the notebook for a different task**
 - New dataset: INSE1920-challenge-ASSIGNMENT.csv
 - **Create new train and test partitions**
 - **Try a different Weka classifier:** Random Forest, NB, SVM, NN, ...
 - Add a paragraph (text cell) at the end briefly **summarizing/discussing the obtained results**
 - (captions are allowed)
 - **Save the .ipynb notebook and submit it to**
fernando.fernandezm@upm.es
with the following subject:
INSE-SmartHealth-Team_XX-Pepito_y_Paquito




New dataset: classes


Running
in place




Step up with
one knee raise



Squat



Hip circles



Alternative assignment (OPTIONAL)



- **Design your own task** by defining a new inventory of gestures:
 - Should be **repetitive movements** (to simplify the recording)
 - For example, “swimming” skills:
 - Backstrokes, Front crawl, Breaststroke, Butterfly, ...
- **Record your own dataset** at home using the “Suite” app
 - Android and iOS versions available
 - Returns a .csv for every recording session
 - Additional libraries available (sampling rate and windowing)
 - At least 1 minute per class, both members
 - Same device, same hand
- **Use the notebook to train and evaluate the system**





Get raw data



- **Physics Toolbox Suite:** Mobile application for recording raw data from smartphone sensors.
- It also let send/upload information to different platforms like Drive, Whatsapp...
- Sampling frequency: 200 Hz

Files
managment

Start/Stop data recording

Plot customization options



Contents

Introduction
HW: Smartphone
SW architecture
Assignment


Weka Experimenter



Weka Experimenter



1. Start Weka



Click the “Experimenter” button to launch the Weka Experimenter.


The Weka Experimenter allows you to design your own experiments of running algorithms on datasets, run the experiments and analyze the results. It’s a powerful tool.



Weka Experimenter



2. Design Experiment



Click the “New” button to create a new experiment configuration.

Test Options


The experimenter configures the test options for you with sensible defaults. The experiment is configured to use Cross Validation with 10 folds. It is a “Classification” type problem and each algorithm + dataset combination is run 10 times (iteration control).



Weka Experimenter



2. Design Experiment



In the “Datasets” select click the “Add new...” button.

Open the “data” directory and choose the “aux_features_raspi.arff” dataset (“iris.arff” in the figure).

Now let’s choose 3 algorithms to run our dataset.



Weka Experimenter



- **ZeroR** is the simplest algorithm we can run. It picks the class value that is the majority in the dataset and gives that for all predictions. It is good to have ZeroR as a baseline that we demand algorithms to outperform.
- **OneR** is like our second simplest algorithm. It picks one attribute that best correlates with the class value and splits it up to get the best prediction accuracy it can.
- **J48** is decision tree algorithm. It is an implementation of the C4.8 algorithm in Java ("J" for Java and 48 for C4.8).

Click "Add new..." in the "Algorithms" section.
Click the "Choose" button.
Click "ZeroR" under the "rules" selection.

Click "Add new..." in the "Algorithms" section.
Click the "Choose" button.
Click "OneR" under the "rules" selection.

Click "Add new..." in the "Algorithms" section.
Click the "Choose" button.
Click "J48" under the "trees" selection.



Weka Experimenter



2. Design Experiment

The screenshot shows the 'Weka Experiment Environment' window with the 'Setup' tab selected. The 'Experiment Configuration Mode' is set to 'Simple'. The 'Results Destination' is 'ARFF file' with a 'Filename' field and a 'Browse...' button. The 'Experiment Type' is 'Cross-validation' with 'Number of folds' set to '10'. The 'Classification' radio button is selected. The 'Iteration Control' shows 'Number of repetitions' set to '10' and 'Data sets first' selected. The 'Datasets' section has 'Add new...', 'Edit selected...', and 'Delete selected...' buttons, with a list containing '/Applications/weka-3-6-10/data/iris.arff'. The 'Algorithms' section has 'Add new...', 'Edit selected...', and 'Delete selected' buttons, with a list containing 'ZeroR', 'OneR -B 6', and 'J48 -C 0.25 -M 2'. At the bottom, there are 'Load options...', 'Save options...', 'Up', and 'Down' buttons, and a 'Notes' field.


We are ready to run our experiment!!!!



Weka Experimenter



3. Run Experiment



Click the “Run” tab at the top of the screen.

This tab is the control panel for running the currently configured experiment.


Click the big “Start” button to start the experiment and watch the “Log” and “Status” sections to keep an eye on how it is doing.

Given that the dataset is small and the algorithms are fast, the experiment should complete in seconds.

Weka Experimenter



4. Review Results



Click the “Analyse” tab at the top of the screen.

This will open up the experiment results analysis panel.


Click the “Experiment” button in the “Source” section to load the results from the current experiment.



Weka Experimenter



4. Review Results: Algorithm Rank



The first thing we want to know is which algorithm was the best. We can do that by ranking the algorithms by the number of times a given algorithm beat the other algorithms.

Click the “Select” button for the “Test base” and choose “Ranking”.


Now Click the “Perform test” button.



Weka Experimenter



4. Review Results: Algorithm Rank



The ranking table shows the number of statistically significant wins each algorithm has had against all other algorithms on the dataset. A win, means an accuracy that is better than the accuracy of another algorithm and that the difference was statistically significant.


We can see that both J48 and OneR have one win each and that ZeroR has two losses. This is good, it means that OneR and J48 are both potentially contenders outperforming out baseline of ZeroR.



Weka Experimenter



4. Review Results: Algorithm Accuracy



Next we want to know what scores the algorithms achieved.

Click the “Select” button for the “Test base” and choose the “ZeroR” algorithm in the list and click the “Select” button.

Click the check-box next to “Show std. deviations”.

Now click the “Perform test” button.



Weka Experimenter



4. Review Results: Algorithm Accuracy

The screenshot shows the Weka Experiment Environment window. The 'Source' tab is active, displaying 'Got 300 results'. The 'Configure test' section on the left includes settings for 'Testing with' (Paired T-Test), 'Row' (Select), 'Column' (Select), 'Comparison field' (Percent_correct), 'Significance' (0.05), 'Sorting (asc.) by' (<default>), 'Test base' (Select), 'Displayed Columns' (Select), 'Show std. deviations' (checked), and 'Output Format' (Select). The 'Test output' section on the right displays the following information:

```
Tester: weka.experiment.PairedCorrectedTTester
Analysing: Percent_correct
Datasets: 1
Resultsets: 3
Confidence: 0.05 (two tailed)
Sorted by: -
Date: 2/19/14 5:42 AM
```

Dataset	(1) rules.ZeroR	(2) rules.OneR	(3) trees.J48
iris	(100) 33.33(0.00)	92.53(5.47) v	94.73(5.30) v
	(v/ /*)	(1/0/0)	(1/0/0)

Key:

```
(1) rules.ZeroR '' 48055541465867954
(2) rules.OneR '-B 6' -3459427003147861443
(3) trees.J48 '-C 0.25 -M 2' -217733168393644444
```

The 'Result list' at the bottom shows three entries:

- 05:37:13 - Percent_correct - Ranking
- 05:37:40 - Percent_correct - trees.J48 -C 0.25
- 05:42:44 - Percent_correct - rules.ZeroR "480

In the “Test output” we can see a table with the results for 3 algorithms. Each algorithm was run 10 times on the dataset and the accuracy reported is the mean and the standard deviation in brackets of those 10 runs.

We can see that both the OneR and J48 algorithms have a little “v” next to their results. This means that the difference in the accuracy for these algorithms compared to ZeroR is statistically significant.

We can also see that the accuracy for these algorithms compared to ZeroR is high, so we can say that these two algorithms achieved a statistically significantly better result than the ZeroR baseline.



References

- Python to GNU Octave bridge
<https://pypi.python.org/pypi/oct2py>
- Python wrapper for the Weka Machine Learning Workbench
<https://pypi.python.org/pypi/python-weka-wrapper>
- WEKA in the Ecosystem for Scientific Computing
http://www.cs.waikato.ac.nz/~eibe/WEKA_Ecosystem.ppt
- Introduction to Google Colab:
https://colab.research.google.com/notebooks/intro.ipynb#scrollTo=5fCEDCU_qrC0

Questions?



